Reichert Clinic-Chek Digital Refractometer



Powerful multi-scale performance in the palm of your hand!

Compact, Powerful, Dependable, Accurate

The Reichert Clinic-Chek is a multi-scale clinical refractometer ideally suited for the medical clinic, laboratory, and for athletic/wrestling hydration measurement. Three measurement scales in one model allows more flexibility and cost savings to you. The Clinic-Chek is more accurate than hydrometers and test strips.

Measures:

- Urine Specific Gravity (USG) from 1.000 to 1.040
- Serum Protein Concentration 1.6 to 15.0g/100mL
- Refractive Index from 1.3330 to 1.4465.

Clinic-Chek raises the bar for performance and value.

The Clinic-Chek is easy to operate and provides quick results. Its LCD digital display is easy to read and removes the subjectivity of shadow line interpretation across a small scale used in conventional handheld refractometers. Reading errors are eliminated with this easy-to-use clinical refractometer.

Measure body fluids (Urine and Blood Serum) without placing the sample near your face as required with other optical manual models. This measurement is safer and more convenient.

Reichert Clinic-Chek Technical Specifications		
Catalog Number	13940021	
USG Range/Accuracy	1.000 to 1.040/0.005	
Serum Protein Range/Accuracy	1.6 to 15.0/0.1	
Refractive Index Range/Accuracy	1.3330 to 1.4465 / 0.0002	
Calibration	Water	
Prism Material	Glass	
Dimensions (WxDxL)	2.1 x 1.1 x 3.9 inches (54 x 27 x 100 mm)	
Weight	3.5 ounces (100 grams)	
Power	10,000 measurements/2xAAA	
Warranty	1 Year	
Ratings & Compliance	IP65 Dustproof/Water Resistant, CE, RoHS, WEEE compliant	



Save Money, Time, and Increase Productivity!

- Clinic-Chek features a 5-digit, full view LCD screen that displays all pertinent measurement data in one view. This feature reduces the costs associated with reading errors and eliminates scrolling values.
 Since the instrument is programmed with multiple measurement scales, only one unit is needed.
- Clinic-Chek has an affordable cost of ownership.
 You can read over 10,000 measurements on the two
 AAA batteries, which are included. The instrument
 also has an automatic shutdown feature that preserves
 battery life.



Analytical Instruments • Refractometers

Innovative precision instruments for over a century

Reichert Provides Real World Solutions and Outstanding Value to Every Customer



Part #	Description	
13107000	AR700 AUTOMATIC REFRACTOMETER	
13107700	AR70 AUTOMATIC REFRACTOMETER	
13106000	AR600 AUTOMATIC REFRACTOMETER	
13106600	AR60 AUTOMATIC REFRACTOMETER	
13980000	r ² i300 AUTOMATIC REFRACTOMETER	



Part #	Description
1310499SA	ARIAS OPTIMATRIX 500 SEMI-AUTOMATIC
	ABBE REFRACTOMETER
1310488M	Abbe MARK III REFRACTOMETER









Analytical Instruments • Refractometers

Corporate Office

Reichert Technologies 3362 Walden Avenue Buffalo, New York 14043 USA Tel. +1 716-686-4500 Fax. +1 716-686-4545 Toll Free USA 1-888-849-8955

e-mail: reichertai.refractometers@ametek.com www.reichertai.com

European Service Center

Carl-von-Linde-Str. 42 85716 Unterschleissheim / Munich Germany Tel: +49 (89) 315 8911 0 Fax: +49 (89) 315 891 99

Reichert Technologies is a division of

